



## Stance Analysis – A Tool for Early Lameness Detection

### ■ What is Stance Analysis?

Stance analysis is a diagnostic procedure that shows how an animal's weight is distributed across each of their limbs. By rapidly capturing multiple snapshots of weight bearing in each limb, and comparing the readings to standard weight-bearing data, the clinician can quickly see if a patient is exhibiting lameness and isolate any issue to specific limb/limbs.



A tool called a “Stance Analyzer” can be used to quickly and easily perform the stance analysis. To perform a stance analysis, the operator places the patient on the special mat, with one paw in each quadrant, and then uses a remote to capture multiple weight bearing measurements that are averaged for a final reading. The results can then be saved and reviewed as part of a wellness, and or lameness assessment.



## ■ What Are the Benefits of Stance Analysis?

Stance analysis offers several benefits, for both rehabilitation-focused and general veterinary practice.

### 1. Earlier Lameness Detection

Musculoskeletal conditions are often difficult to detect until an animal starts showing outward signs, such as limping. By using a Stance Analyzer, veterinary professionals can objectively see when there's a problem before it manifests as modified behavior. It can also help pinpoint the location of the problem (i.e. left rear limb), enabling faster diagnosis, and development of tailored treatment plans. Additionally, a Stance Analyzer provides early lameness detection comparable to a pressure sensitive walkway<sup>1</sup>, but has a smaller footprint and lower capital outlay.


### 2. Better Client Communication

While some pet owners take their veterinarian's treatment recommendations without hesitation, many need a bit more convincing – especially when they haven't noticed anything strange in their pet's behavior. Having objective data allows the clinician to show the owner that a problem exists. Their pet may be in pain, and early intervention will save future suffering and expense. Using objective data also promotes a willingness to begin a recommended treatment regimen right away, and drives greater compliance.

### 3. Clearer Progress Tracking

For patients who present unable to walk or with a pronounced limp, improvement may be easy for both the veterinarian and the pet owner to see as treatment progresses. For less subtle changes, however, objective data is critical to ongoing compliance. By repeating stance analysis, and comparing the results over time, it is easy to see improvement as the weight-bearing percentages normalize.

1. Clough, William T. et al. "Sensitivity and Specificity of a Weight Distribution Platform for the Detection of Objective Lameness and Orthopaedic Disease." *Veterinary and Comparative Orthopaedics and Traumatology*, Vol 252, No. 6, 2018, pp. 701-709.



## ■ When Should Stance Analysis Be Used?

For maximum benefit, stance analysis should be performed during wellness exams, problem appointments/diagnostic visits, and treatment re-checks.

Using the Stance Analyzer as a part of a standard wellness exam allows veterinary professionals to identify lameness that may not yet be obvious. It can also be used as a scale for routine weight measurement.

When used during a problem appointment, stance analysis helps to pinpoint where a problem is originating. It additionally provides objective pain/lameness data that may be viewed, stored, or printed and given to the pet owner.



During the course of treatment, stance analysis should be used to measure progress. Pet owners appreciate being able to see objective data of progress, both for the relief it provides to see their pet is improving, but also to justify the investment they've made in the specific treatment.

## The Companion Stance Analyzer 3 Daily Applications



### Wellness Exams



#### SCALE

Replace your scale with the Stance Analyzer for weight measurement at every visit



#### TRACK TRENDS

Take weight bearing data at every visit to identify weight bearing trends



#### PREVENTION

Identify abnormalities earlier, recommend additional plans of care



### Diagnostic Visits



#### DIAGNOSE

Use the data from the Stance Analyzer to establish a base line and to determine where to look



#### SERVICES

The data validates which services you are prescribing and increases client compliance



#### POINT-OF-CARE

An instant results diagnostic tool that the pet owner can participate in



### Treatment Re-Checks



#### QUANTIFY

Data tells you when a pet has improved by an objective % and is finished with their treatment



#### VALIDATE

Objective data validates the success of your treatment plan and helps with compliance



#### PIVOT

When a treatment is not producing the desired outcomes use the data to support plan changes

For more information or to book a demonstration of a Companion Stance Analyzer, please contact Georgina MacPhail at 01646 603878 or [georginam@litecure.com](mailto:georginam@litecure.com).



[www.companionanimalhealth.com](http://www.companionanimalhealth.com)